

# The Choice for Electromedical Procedures

## Spectra® 360 Electrode Gel: Safe - Salt-free - Bacteriostatic

Spectra 360's totally salt-free/chloride-free electrode gel is non-corrosive to electrodes and improves the conductive properties of rubber electrodes.

Its exclusive STAY-WET® formula is optimum for prolonged use in ECG, TENS, pediatric, and all long term procedures—requiring no re-application.

### Spectra 360 Electrode Gel...

- Bacteriostatic formula minimizes risk of infection
- Reduces skin resistance
- Non-irritating
- Does not stain clothing or linens
- Leaves no caking or residue



ISO 13485:2003



Parker Laboratories, Inc.

The sound choice in patient care.™

800.631.8888

[www.parkerlabs.com](http://www.parkerlabs.com)

## SPECTRA® 360 ELECTRODE GEL

The only salt-free and chloride-free electrically conductive gel, recommended for all electromedical procedures, except defibrillation. Salt-free characteristics make it particularly suitable for long-term applications. Spectra 360 differs significantly from all other electrically conductive media... it works by wetting the skin, thereby reducing skin resistance.

### FEATURES

- Non-gritty STAY-WET® formula allows for prolonged use without re-application
- Product of choice for ECG, TENS, monitoring and pediatrics
- Non-irritating, hypoallergenic, bacteriostatic
- Salt-free, no sodium ion transfer
- Product of choice for conductive rubber/carbon electrodes
- Not recommended for defibrillation

Product #	Package Size
12-02	60 g (2 oz.) tube, 12 per box
12-08	250 g (8.5 oz.) tube, 12 per box



ISO 13485:2003



Parker Laboratories, Inc.

The sound choice in patient care.™

800.631.8888

www.parkerlabs.com